

Get the best from your Sonicare DiamondClean - Quick Guide

Charging your DiamondClean

Your Sonicare DiamondClean comes with two charging methods, offering you charging options at home and away. In your box you will find the following:



Using the charger base

Simply place the handle in the glass and place the glass on the charger base. Plug the charger into a 2-pin shaver socket, or a standard socket using a shaver plug adapter.

Using the travel charger

While travelling you can use the charging travel case in two ways.

1. To charge via wall outlet:

Place the handle in the charging travel case, plug one end of the USB cord into the travel case, and plug the USB end into the USB wall adapter and plug this into an electrical (shaver) 2-pin plug socket. The charge indicator on the handle flashes green until fully charged.

2 To charge via USB:

Insert the USB end of the charger cable into any USB port on your computer.

Use only the Philips USB wall adapter and USB cable provided.

About the 2-pin plug

Your Sonicare DiamondClean is supplied with a UK 'shaver' plug which is mandatory for bathroom appliances. This plug should only be plugged into a shaver plug socket, which you may find in your bathroom. If you do not have a shaver plug socket in your bathroom you can buy an adaptor specially designed for shaver plugs, to convert and use in a standard UK 3-pin plug socket.

Getting the most out of your battery

Your Sonicare DiamondClean has been designed to last for 3 weeks on a single full charge, based on two minutes of brushing, twice a day. The actual number of brushings you get will depend on how often you brush, the number of people using the toothbrush and the length of time the handle has been charging. Ensure you charge your DiamondClean for a full 24 hours before first use. The charger flashes green until the handle is fully charged.

Your DiamondClean brush head needs space

You may notice a small gap between the top of the toothbrush handle and the base of the brush head. This gap is necessary, as the brush head must be free to move in order to deliver the correct brushing action needed to clean the teeth.

Cleaning and storing your DiamondClean

Do not leave your DiamondClean in water – wipe the handle dry after each use and clean it periodically. To clean, rinse the handle and the brush head. Wipe down all surfaces with a soft cloth and mild soap, if desired. Do not use sharp objects to clean.

If you are not going to use your DiamondClean for an extended period of time, unplug the charging base, clean all parts and store them in a cool, dry place away from direct sunlight.